



MECH ORTHO PVT. LTD.

Reconstruct foundation of human body

BONE SCREW



online shop

mechorthoebuy.com

Mobile App



Online Store




Website




www.mechortho.com

1.50mm CORTICAL SCREW




THREAD DIA.	CORE DIA.	HEAD DIA	HEXAGONAL SOCKET	TPI	
1.50mm	1.0mm	3.0mm	1.50mm	40	
TOTAL LENGTH	CATALOGUE NO.		TOTAL LENGTH	CATALOGUE NO.	
	SS.316L	TITANIUM		SS.316L	TITANIUM
6mm	1.001.006	Ti.1.001.006	20mm	1.001.020	Ti.1.001.020
8mm	1.001.008	Ti.1.001.008	22mm	1.001.022	Ti.1.001.022
10mm	1.001.010	Ti.1.001.010	24mm	1.001.024	Ti.1.001.024
12mm	1.001.012	Ti.1.001.012	26mm	1.001.026	Ti.1.001.026
14mm	1.001.014	Ti.1.001.014	28mm	1.001.028	Ti.1.001.028
16mm	1.001.016	Ti.1.001.016	30mm	1.001.030	Ti.1.001.030
18mm	1.001.018	Ti.1.001.018			

2.70mm CORTICAL SCREW




THREAD DIA.	CORE DIA.	HEAD DIA	HEXAGONAL SOCKET	TPI	
2.70mm	1.90mm	5.0mm	2.50mm	24	
TOTAL LENGTH	CATALOGUE NO.		TOTAL LENGTH	CATALOGUE NO.	
	SS.316L	TITANIUM		SS.316L	TITANIUM
6mm	1.002.006	Ti.1.002.006	24mm	1.002.024	Ti.1.002.024
8mm	1.002.008	Ti.1.002.008	26mm	1.002.026	Ti.1.002.026
10mm	1.002.010	Ti.1.002.010	28mm	1.002.028	Ti.1.002.028
12mm	1.002.012	Ti.1.002.012	30mm	1.002.030	Ti.1.002.030
14mm	1.002.014	Ti.1.002.014	32mm	1.002.032	Ti.1.002.032
16mm	1.002.016	Ti.1.002.016	34mm	1.002.034	Ti.1.002.034
18mm	1.002.018	Ti.1.002.018	36mm	1.002.036	Ti.1.002.036
20mm	1.002.020	Ti.1.002.020	38mm	1.002.038	Ti.1.002.038
22mm	1.002.022	Ti.1.002.022	40mm	1.002.040	Ti.1.002.040

3.50mm CORTICAL SCREW




THREAD DIA.	CORE DIA.	HEAD DIA.	HEXAGONAL SOCKET	TPI	
3.50mm	2.40mm	6.0mm	2.50mm	20	
TOTAL LENGTH	CATALOGUE NO.		TOTAL LENGTH	CATALOGUE NO.	
	SS.316L	TITANIUM		SS.316L	TITANIUM
10mm	1.003.010	Ti.1.003.010	32mm	1.003.032	Ti.1.003.032
12mm	1.003.012	Ti.1.003.012	34mm	1.003.034	Ti.1.003.034
14mm	1.003.014	Ti.1.003.014	36mm	1.003.036	Ti.1.003.036
16mm	1.003.016	Ti.1.003.016	38mm	1.003.038	Ti.1.003.038
18mm	1.003.018	Ti.1.003.018	40mm	1.003.040	Ti.1.003.040
20mm	1.003.020	Ti.1.003.020	42mm	1.003.042	Ti.1.003.042
22mm	1.003.022	Ti.1.003.022	45mm	1.003.045	Ti.1.003.045
24mm	1.003.024	Ti.1.003.024	50mm	1.003.050	Ti.1.003.050
26mm	1.003.026	Ti.1.003.026	55mm	1.003.055	Ti.1.003.055
28mm	1.003.028	Ti.1.003.028	60mm	1.003.060	Ti.1.003.060
30mm	1.003.030	Ti.1.003.030			

4.50mm CORTICAL SCREW




THREAD DIA.	CORE DIA.	HEAD DIA.	HEXAGONAL SOCKET	TPI	
4.50mm	3.0mm	8.0mm	3.50mm	14	
TOTAL LENGTH	CATALOGUE NO.		TOTAL LENGTH	CATALOGUE NO.	
	SS.316L	TITANIUM		SS.316L	TITANIUM
14mm	1.004.014	Ti.1.004.014	46mm	1.004.046	Ti.1.004.046
16mm	1.004.016	Ti.1.004.016	48mm	1.004.048	Ti.1.004.048
18mm	1.004.018	Ti.1.004.018	50mm	1.004.050	Ti.1.004.050
20mm	1.004.020	Ti.1.004.020	52mm	1.004.052	Ti.1.004.052
22mm	1.004.022	Ti.1.004.022	54mm	1.004.054	Ti.1.004.054
24mm	1.004.024	Ti.1.004.024	56mm	1.004.056	Ti.1.004.056
26mm	1.004.026	Ti.1.004.026	58mm	1.004.058	Ti.1.004.058
28mm	1.004.028	Ti.1.004.028	60mm	1.004.060	Ti.1.004.060
30mm	1.004.030	Ti.1.004.030	62mm	1.004.062	Ti.1.004.062
32mm	1.004.032	Ti.1.004.032	64mm	1.004.064	Ti.1.004.064
34mm	1.004.034	Ti.1.004.034	66mm	1.004.066	Ti.1.004.066
36mm	1.004.036	Ti.1.004.036	68mm	1.004.068	Ti.1.004.068
38mm	1.004.038	Ti.1.004.038	70mm	1.004.070	Ti.1.004.070
40mm	1.004.040	Ti.1.004.040	72mm	1.004.072	Ti.1.004.072
42mm	1.004.042	Ti.1.004.042	74mm	1.004.074	Ti.1.004.074
44mm	1.004.044	Ti.1.004.044			

2.70mm CORTICAL SCREW SELF TAPPING




THREAD DIA.	CORE DIA.	HEAD DIA.	HEXAGONAL SOCKET	TPI	
2.70mm	1.90mm	5.0mm	2.50mm	24	
TOTAL LENGTH	CATALOGUE NO.		TOTAL LENGTH	CATALOGUE NO.	
	SS.316L	TITANIUM		SS.316L	TITANIUM
6mm	1.005.006	Ti.1.005.006	24mm	1.005.024	Ti.1.005.024
8mm	1.005.008	Ti.1.005.008	26mm	1.005.026	Ti.1.005.026
10mm	1.005.010	Ti.1.005.010	28mm	1.005.028	Ti.1.005.028
12mm	1.005.012	Ti.1.005.012	30mm	1.005.030	Ti.1.005.030
14mm	1.005.014	Ti.1.005.014	32mm	1.005.032	Ti.1.005.032
16mm	1.005.016	Ti.1.005.016	34mm	1.005.034	Ti.1.005.034
18mm	1.005.018	Ti.1.005.018	36mm	1.005.036	Ti.1.005.036
20mm	1.005.020	Ti.1.005.020	38mm	1.005.038	Ti.1.005.038
22mm	1.005.022	Ti.1.005.022	40mm	1.005.040	Ti.1.005.040

3.50mm CORTICAL SCREW SELF TAPPING




THREAD DIA.	CORE DIA.	HEAD DIA.	HEXAGONAL SOCKET	TPI	
3.50mm	2.40mm	6.0mm	2.50mm	20	
TOTAL LENGTH	CATALOGUE NO.		TOTAL LENGTH	CATALOGUE NO.	
	SS.316L	TITANIUM		SS.316L	TITANIUM
10mm	1.006.010	Ti.1.006.010	32mm	1.006.032	Ti.1.006.032
12mm	1.006.012	Ti.1.006.012	34mm	1.006.034	Ti.1.006.034
14mm	1.006.014	Ti.1.006.014	36mm	1.006.036	Ti.1.006.036
16mm	1.006.016	Ti.1.006.016	38mm	1.006.038	Ti.1.006.038
18mm	1.006.018	Ti.1.006.018	40mm	1.006.040	Ti.1.006.040
20mm	1.006.020	Ti.1.006.020	42mm	1.006.042	Ti.1.006.042
22mm	1.006.022	Ti.1.006.022	45mm	1.006.044	Ti.1.006.044
24mm	1.006.024	Ti.1.006.024	50mm	1.006.046	Ti.1.006.046
26mm	1.006.026	Ti.1.006.026	55mm	1.006.048	Ti.1.006.048
28mm	1.006.028	Ti.1.006.028	60mm	1.006.050	Ti.1.006.050
30mm	1.006.030	Ti.1.006.030			

4.50mm CORTICAL SCREW SELF TAPPING




THREAD DIA.	CORE DIA.	HEAD DIA.	HEXAGONAL SOCKET	TPI	
4.50mm	3.0mm	8.0mm	3.50mm	14	
TOTAL LENGTH	CATALOGUE NO.		TOTAL LENGTH	CATALOGUE NO.	
	SS.316L	TITANIUM		SS.316L	TITANIUM
14mm	1.007.014	Ti.1.007.014	46mm	1.007.046	Ti.1.007.046
16mm	1.007.016	Ti.1.007.016	48mm	1.007.048	Ti.1.007.048
18mm	1.007.018	Ti.1.007.018	50mm	1.007.050	Ti.1.007.050
20mm	1.007.020	Ti.1.007.020	52mm	1.007.052	Ti.1.007.052
22mm	1.007.022	Ti.1.007.022	54mm	1.007.054	Ti.1.007.054
24mm	1.007.024	Ti.1.007.024	56mm	1.007.056	Ti.1.007.056
26mm	1.007.026	Ti.1.007.026	58mm	1.007.058	Ti.1.007.058
28mm	1.007.028	Ti.1.007.028	60mm	1.007.060	Ti.1.007.060
30mm	1.007.030	Ti.1.007.030	62mm	1.007.062	Ti.1.007.062
32mm	1.007.032	Ti.1.007.032	64mm	1.007.064	Ti.1.007.064
34mm	1.007.034	Ti.1.007.034	66mm	1.007.066	Ti.1.007.066
36mm	1.007.036	Ti.1.007.036	68mm	1.007.068	Ti.1.007.068
38mm	1.007.038	Ti.1.007.038	70mm	1.007.070	Ti.1.007.070
40mm	1.007.040	Ti.1.007.040	72mm	1.007.072	Ti.1.007.072
42mm	1.007.042	Ti.1.007.042	74mm	1.007.074	Ti.1.007.074
44mm	1.007.044	Ti.1.007.044			

2.70mm LOCKING HEAD SELF TAPPING SCREW




THREAD DIA.	CORE DIA.	HEAD DIA.	HEXAGONAL SOCKET	TPI	
2.70mm	1.90mm	3.50mm	2.0mm	48	
TOTAL LENGTH	CATALOGUE NO.		TOTAL LENGTH	CATALOGUE NO.	
	SS.316L	TITANIUM		SS.316L	TITANIUM
10mm	1.008.010	Ti.1.008.010	26mm	1.008.026	Ti.1.008.026
12mm	1.008.012	Ti.1.008.012	28mm	1.008.028	Ti.1.008.028
14mm	1.008.014	Ti.1.008.014	30mm	1.008.030	Ti.1.008.030
16mm	1.008.016	Ti.1.008.016	32mm	1.008.032	Ti.1.008.032
18mm	1.008.018	Ti.1.008.018	34mm	1.008.034	Ti.1.008.034
20mm	1.008.020	Ti.1.008.020	36mm	1.008.036	Ti.1.008.036
22mm	1.008.022	Ti.1.008.022	38mm	1.008.038	Ti.1.008.038
24mm	1.008.024	Ti.1.008.024	40mm	1.008.040	Ti.1.008.040

3.50mm LOCKING HEAD SELF TAPPING SCREW




THREAD DIA.	CORE DIA.	HEAD DIA.	HEXAGONAL SOCKET	TPI	
3.50mm	2.60mm	4.80mm	2.50mm	32	
TOTAL LENGTH	CATALOGUE NO.		TOTAL LENGTH	CATALOGUE NO.	
	SS.316L	TITANIUM		SS.316L	TITANIUM
10mm	1.009.010	Ti.1.009.010	28mm	1.009.028	Ti.1.009.028
12mm	1.009.012	Ti.1.009.012	30mm	1.009.030	Ti.1.009.030
14mm	1.009.014	Ti.1.009.014	35mm	1.009.035	Ti.1.009.035
16mm	1.009.016	Ti.1.009.016	40mm	1.009.040	Ti.1.009.040
18mm	1.009.018	Ti.1.009.018	45mm	1.009.045	Ti.1.009.045
20mm	1.009.020	Ti.1.009.020	50mm	1.009.050	Ti.1.009.050
22mm	1.009.022	Ti.1.009.022	55mm	1.009.055	Ti.1.009.055
24mm	1.009.024	Ti.1.009.024	60mm	1.009.060	Ti.1.009.060
26mm	1.009.026	Ti.1.009.026			

5.0mm LOCKING HEAD SELF TAPPING SCREW




THREAD DIA.	CORE DIA.	HEAD DIA.	HEXAGONAL SOCKET	TPI	
5.0mm	4.30mm	6.80mm	3.50mm	26	
TOTAL LENGTH	CATALOGUE NO.		TOTAL LENGTH	CATALOGUE NO.	
	SS.316L	TITANIUM		SS.316L	TITANIUM
14mm	1.010.014	Ti.1.010.014	42mm	1.010.042	Ti.1.010.042
16mm	1.010.016	Ti.1.010.016	44mm	1.010.044	Ti.1.010.044
18mm	1.010.018	Ti.1.010.018	46mm	1.010.046	Ti.1.010.046
20mm	1.010.020	Ti.1.010.020	48mm	1.010.048	Ti.1.010.048
22mm	1.010.022	Ti.1.010.022	50mm	1.010.050	Ti.1.010.050
24mm	1.010.024	Ti.1.010.024	55mm	1.010.055	Ti.1.010.055
26mm	1.010.024	Ti.1.010.026	60mm	1.010.060	Ti.1.010.060
28mm	1.010.026	Ti.1.010.028	65mm	1.010.065	Ti.1.010.065
30mm	1.010.028	Ti.1.010.030	70mm	1.010.070	Ti.1.010.070
32mm	1.010.030	Ti.1.010.032	75mm	1.010.075	Ti.1.010.075
34mm	1.010.032	Ti.1.010.034	80mm	1.010.080	Ti.1.010.080
36mm	1.010.034	Ti.1.010.036	85mm	1.010.085	Ti.1.010.085
38mm	1.010.038	Ti.1.010.038	90mm	1.010.090	Ti.1.010.090
40mm	1.010.040	Ti.1.010.040			

4.0mm CANCELLOUS SCREW - Short Thread




THREAD DIA.	CORE DIA.	HEAD DIA	HEXAGONAL SOCKET	TPI	
4.0mm	1.90mm	6.0mm	2.50mm	14	
TOTAL LENGTH	CATALOGUE NO.		TOTAL LENGTH	CATALOGUE NO.	
	SS.316L	TITANIUM		SS.316L	TITANIUM
10mm	1.011.010	1.011.010	28mm	1.011.028	1.011.028
12mm	1.011.012	1.011.012	30mm	1.011.030	1.011.030
14mm	1.011.014	1.011.014	35mm	1.011.035	1.011.035
16mm	1.011.016	1.011.016	40mm	1.011.040	1.011.040
18mm	1.011.018	1.011.018	45mm	1.011.045	1.011.045
20mm	1.011.020	1.011.020	50mm	1.011.050	1.011.050
22mm	1.011.022	1.011.022	55mm	1.011.055	1.011.055
24mm	1.011.024	1.011.024	60mm	1.011.060	1.011.060
26mm	1.011.026	1.011.026			

6.50mm Cancellous Screw - 16mm Thread




THREAD DIA.	CORE DIA.	HEAD DIA	HEXAGONAL SOCKET	TPI	
6.50mm	3.0mm	8.0mm	3.50mm	9	
TOTAL LENGTH	CATALOGUE NO.		TOTAL LENGTH	CATALOGUE NO.	
	SS.316L	TITANIUM		SS.316L	TITANIUM
30.0	1.012.030	Ti.1.012.030	75.0	1.012.075	Ti.1.012.075
35.0	1.012.035	Ti.1.012.035	80.0	1.012.080	Ti.1.012.080
40.0	1.012.040	Ti.1.012.040	85.0	1.012.085	Ti.1.012.085
45.0	1.012.045	Ti.1.012.045	90.0	1.012.090	Ti.1.012.090
50.0	1.012.050	Ti.1.012.050	95.0	1.012.095	Ti.1.012.095
55.0	1.012.055	Ti.1.012.055	100.0	1.012.100	Ti.1.012.100
60.0	1.012.060	Ti.1.012.060	105.0	1.012.105	Ti.1.012.105
65.0	1.012.065	Ti.1.012.065	110.0	1.012.110	Ti.1.012.110
70.0	1.012.070	Ti.1.012.070			

6.50mm Cancellous Screw - 16mm Thread




THREAD DIA.	CORE DIA.	HEAD DIA	HEXAGONAL SOCKET	TPI	
6.50mm	3.0mm	8.0mm	3.50mm	9	
TOTAL LENGTH	CATALOGUE NO.		TOTAL LENGTH	CATALOGUE NO.	
	SS.316L	TITANIUM		SS.316L	TITANIUM
40.0	1.013.040	Ti.1.013.040	80.0	1.012.080	Ti.1.012.080
45.0	1.013.045	Ti.1.013.045	85.0	1.012.085	Ti.1.012.085
50.0	1.013.050	Ti.1.013.050	90.0	1.012.090	Ti.1.012.090
55.0	1.013.055	Ti.1.013.055	95.0	1.012.095	Ti.1.012.095
60.0	1.013.060	Ti.1.013.060	100.0	1.012.100	Ti.1.012.100
65.0	1.013.065	Ti.1.013.065	105.0	1.012.105	Ti.1.012.105
70.0	1.012.070	Ti.1.013.070	110.0	1.012.110	Ti.1.012.110
75.0	1.012.075	Ti.1.013.075			

4.0mm CANCELLOUS SCREW - Full thread




THREAD DIA.	CORE DIA.	HEAD DIA	HEXAGONAL SOCKET	TPI	
4.0mm	1.90mm	6.0mm	2.50mm	14	
TOTAL LENGTH	CATALOGUE NO.		TOTAL LENGTH	CATALOGUE NO.	
	SS.316L	TITANIUM		SS.316L	TITANIUM
10mm	1.014.010	Ti.1.014.010	28mm	1.014.028	Ti.1.014.028
12mm	1.014.012	Ti.1.014.012	30mm	1.014.030	Ti.1.014.030
14mm	1.014.014	Ti.1.014.014	35mm	1.014.035	Ti.1.014.035
16mm	1.014.016	Ti.1.014.016	40mm	1.014.040	Ti.1.014.040
18mm	1.014.018	Ti.1.014.018	45mm	1.014.045	Ti.1.014.045
20mm	1.014.020	Ti.1.014.020	50mm	1.014.050	Ti.1.014.050
22mm	1.014.022	Ti.1.014.022	55mm	1.014.055	Ti.1.014.055
24mm	1.014.024	Ti.1.014.024	60mm	1.014.060	Ti.1.014.060
26mm	1.014.026	Ti.1.014.026			

6.50mm Cancellous Screw - Full Thread




THREAD DIA.	CORE DIA.	HEAD DIA.	HEXAGONAL SOCKET	TPI	
6.50mm	3.0mm	8.0mm	3.50mm	9	
TOTAL LENGTH	CATALOGUE NO.		TOTAL LENGTH	CATALOGUE NO.	
	SS.316L	TITANIUM		SS.316L	TITANIUM
25mm	1.015.025	Ti.1.015.025	70mm	1.015.070	Ti.1.015.070
30mm	1.015.030	Ti.1.015.030	75mm	1.015.075	Ti.1.015.075
35mm	1.015.035	Ti.1.015.035	80mm	1.015.080	Ti.1.015.080
40mm	1.015.040	Ti.1.015.040	85mm	1.015.085	Ti.1.015.085
45mm	1.015.045	Ti.1.015.045	90mm	1.015.090	Ti.1.015.090
50mm	1.015.050	Ti.1.015.050	95mm	1.015.095	Ti.1.015.095
55mm	1.015.055	Ti.1.015.055	100mm	1.015.100	Ti.1.015.100
60mm	1.015.060	Ti.1.015.060	105mm	1.015.105	Ti.1.015.105
65mm	1.015.065	Ti.1.015.065	110mm	1.015.110	Ti.1.015.110

4.0mm Cannulated Cancellous Screw - Full Thread




THREAD DIA.	CORE DIA.	HEAD DIA.	HEXAGONAL SOCKET	CANNULATION	
4.0mm	2.60mm	6.0mm	2.50mm	1.40mm	
TOTAL LENGTH	CATALOGUE NO.		TOTAL LENGTH	CATALOGUE NO.	
	SS.316L	TITANIUM		SS.316L	TITANIUM
10mm	1.016.010	Ti.1.016.010	34mm	1.016.034	Ti.1.016.034
12mm	1.016.012	Ti.1.016.012	36mm	1.016.036	Ti.1.016.036
14mm	1.016.014	Ti.1.016.014	38mm	1.016.038	Ti.1.016.038
16mm	1.016.016	Ti.1.016.016	40mm	1.016.040	Ti.1.016.040
18mm	1.016.018	Ti.1.016.018	42mm	1.016.042	Ti.1.016.042
20mm	1.016.020	Ti.1.016.020	44mm	1.016.044	Ti.1.016.044
22mm	1.016.022	Ti.1.016.022	46mm	1.016.046	Ti.1.016.046
24mm	1.016.024	Ti.1.016.024	48mm	1.016.048	Ti.1.016.048
26mm	1.016.026	Ti.1.016.026	50mm	1.016.050	Ti.1.016.050
28mm	1.016.028	Ti.1.016.028	55mm	1.016.055	Ti.1.016.055
30mm	1.016.030	Ti.1.016.030	60mm	1.016.060	Ti.1.016.060
32mm	1.016.032	Ti.1.016.032			

4.0mm Cannulated Cancellous Screw - Short Thread




THREAD DIA.	CORE DIA.	HEAD DIA	HEXAGONAL SOCKET	CANNULATION	
4.0mm	2.60mm	6.0mm	2.50mm	1.40mm	
TOTAL LENGTH	CATALOGUE NO.		TOTAL LENGTH	CATALOGUE NO.	
	SS.316L	TITANIUM		SS.316L	TITANIUM
10mm	1.017.010	Ti.1.017.010	34mm	1.017.034	Ti.1.017.034
12mm	1.017.012	Ti.1.017.012	36mm	1.017.036	Ti.1.017.036
14mm	1.017.014	Ti.1.017.014	38mm	1.017.038	Ti.1.017.038
16mm	1.017.016	Ti.1.017.016	40mm	1.017.040	Ti.1.017.040
18mm	1.017.018	Ti.1.017.018	42mm	1.017.042	Ti.1.017.042
20mm	1.017.020	Ti.1.017.020	44mm	1.017.044	Ti.1.017.044
22mm	1.017.022	Ti.1.017.022	46mm	1.017.046	Ti.1.017.046
24mm	1.017.024	Ti.1.017.024	48mm	1.017.048	Ti.1.017.048
26mm	1.017.026	Ti.1.017.026	50mm	1.017.050	Ti.1.017.050
28mm	1.017.028	Ti.1.017.028	55mm	1.017.055	Ti.1.017.055
30mm	1.017.030	Ti.1.017.030	60mm	1.017.060	Ti.1.017.060
32mm	1.017.032	Ti.1.017.032			

7.0mm Cannulated Cancellous Screw - 16mm Thread



THREAD DIA.	CORE DIA.	HEAD DIA	HEXAGONAL SOCKET	CANNULATION	
7.0mm	4.50mm	8.0mm	3.50mm	2.20mm	
TOTAL LENGTH	CATALOGUE NO.		TOTAL LENGTH	CATALOGUE NO.	
	SS.316L	TITANIUM		SS.316L	TITANIUM
30mm	1.018.030	Ti.1.018.030	80mm	1.018.080	Ti.1.018.080
35mm	1.018.035	Ti.1.018.035	85mm	1.018.085	Ti.1.018.085
40mm	1.018.040	Ti.1.018.040	90mm	1.018.090	Ti.1.018.090
45mm	1.018.045	Ti.1.018.045	95mm	1.018.095	Ti.1.018.095
50mm	1.018.050	Ti.1.018.050	100mm	1.018.100	Ti.1.018.100
55mm	1.018.055	Ti.1.018.055	105mm	1.018.105	Ti.1.018.105
60mm	1.018.060	Ti.1.018.060	110mm	1.018.110	Ti.1.018.110
65mm	1.018.065	Ti.1.018.065	115mm	1.018.115	Ti.1.018.115
70mm	1.018.070	Ti.1.018.070	120mm	1.018.120	Ti.1.018.120
75mm	1.018.075	Ti.1.018.075			

7.0mm Cannulated Cancellous Screw - 32mm Thread



THREAD DIA.	CORE DIA.	HEAD DIA	HEXAGONAL SOCKET	CANNULATION	
7.0mm	4.50mm	8.0mm	3.50mm	2.20mm	
TOTAL LENGTH	CATALOGUE NO.		TOTAL LENGTH	CATALOGUE NO.	
	SS316L	TITANIUM		SS316L	TITANIUM
30mm	1.019.030	Ti.1.019.030	80mm	1.019.080	Ti.1.019.080
35mm	1.019.035	Ti.1.019.035	85mm	1.019.085	Ti.1.019.085
40mm	1.019.040	Ti.1.019.040	90mm	1.019.090	Ti.1.019.090
45mm	1.019.045	Ti.1.019.045	95mm	1.019.095	Ti.1.019.095
50mm	1.019.050	Ti.1.019.050	100mm	1.019.100	Ti.1.019.100
55mm	1.019.055	Ti.1.019.055	105mm	1.019.105	Ti.1.019.105
60mm	1.019.060	Ti.1.019.060	110mm	1.019.110	Ti.1.019.110
65mm	1.019.065	Ti.1.019.065	115mm	1.019.115	Ti.1.019.115
70mm	1.019.070	Ti.1.019.070	120mm	1.019.120	Ti.1.019.120
75mm	1.019.075	Ti.1.019.075			

3.40mm Locking Bolt



THREAD DIA.	CORE DIA.	HEAD DIA	HEXAGONAL SOCKET	TPI	
3.40mm	2.80mm	7.90mm	3.50mm	12	
TOTAL LENGTH	CATALOGUE NO.		TOTAL LENGTH	CATALOGUE NO.	
	SS316L	TITANIUM		SS316L	TITANIUM
20mm	1.020.020	Ti.1.020.020	42mm	1.020.042	Ti.1.020.042
22mm	1.020.022	Ti.1.020.022	44mm	1.020.044	Ti.1.020.044
24mm	1.020.024	Ti.1.020.024	46mm	1.020.046	Ti.1.020.046
26mm	1.020.026	Ti.1.020.026	48mm	1.020.048	Ti.1.020.048
28mm	1.020.028	Ti.1.020.028	50mm	1.020.050	Ti.1.020.050
30mm	1.020.030	Ti.1.020.030	52mm	1.020.052	Ti.1.020.052
32mm	1.020.032	Ti.1.020.032	54mm	1.020.054	Ti.1.020.054
34mm	1.020.034	Ti.1.020.034	56mm	1.020.056	Ti.1.020.056
36mm	1.020.036	Ti.1.020.036	58mm	1.020.058	Ti.1.020.058
38mm	1.020.038	Ti.1.020.038	60mm	1.020.060	Ti.1.020.060
40mm	1.020.040	Ti.1.020.040			

3.90mm Locking Bolt



THREAD DIA.	CORE DIA.	HEAD DIA.	HEXAGONAL SOCKET	TPI	
3.90mm	3.40mm	7.90mm	3.50mm	12	
TOTAL LENGTH	CATALOGUE NO.		TOTAL LENGTH	CATALOGUE NO.	
	SS316L	TITANIUM		SS316L	TITANIUM
20mm	1.021.020	Ti.1.021.020	42mm	1.021.042	Ti.1.021.042
22mm	1.021.022	Ti.1.021.022	44mm	1.021.044	Ti.1.021.044
24mm	1.021.024	Ti.1.021.024	46mm	1.021.046	Ti.1.021.046
26mm	1.021.026	Ti.1.021.026	48mm	1.021.048	Ti.1.021.048
28mm	1.021.028	Ti.1.021.028	50mm	1.021.050	Ti.1.021.050
30mm	1.021.030	Ti.1.021.030	52mm	1.021.052	Ti.1.021.052
32mm	1.021.032	Ti.1.021.032	54mm	1.021.054	Ti.1.021.054
34mm	1.021.034	Ti.1.021.034	56mm	1.021.056	Ti.1.021.056
36mm	1.021.036	Ti.1.021.036	58mm	1.021.058	Ti.1.021.058
38mm	1.021.038	Ti.1.021.038	60mm	1.021.060	Ti.1.021.060
40mm	1.021.040	Ti.1.021.040			

4.90mm Locking Bolt



THREAD DIA.	CORE DIA.	HEAD DIA.	HEXAGONAL SOCKET	TPI	
4.90mm	4.30mm	7.90mm	3.50mm	9	
TOTAL LENGTH	CATALOGUE NO.		TOTAL LENGTH	CATALOGUE NO.	
	SS316L	TITANIUM		SS316L	TITANIUM
24mm	1.022.024	Ti.1.022.024	52mm	1.022.052	Ti.1.022.052
26mm	1.022.026	Ti.1.022.026	54mm	1.022.054	Ti.1.022.054
28mm	1.022.028	Ti.1.022.028	56mm	1.022.056	Ti.1.022.056
30mm	1.022.030	Ti.1.022.030	58mm	1.022.058	Ti.1.022.058
32mm	1.022.032	Ti.1.022.032	60mm	1.022.060	Ti.1.022.060
34mm	1.022.034	Ti.1.022.034	65mm	1.022.065	Ti.1.022.065
36mm	1.022.036	Ti.1.022.036	70mm	1.022.070	Ti.1.022.070
38mm	1.022.038	Ti.1.022.038	75mm	1.022.075	Ti.1.022.075
40mm	1.022.040	Ti.1.022.040	80mm	1.022.080	Ti.1.022.080
42mm	1.022.042	Ti.1.022.042	85mm	1.022.085	Ti.1.022.085
44mm	1.022.044	Ti.1.022.044	90mm	1.022.090	Ti.1.022.090
46mm	1.022.046	Ti.1.022.046	95mm	1.022.095	Ti.1.022.095
48mm	1.022.048	Ti.1.022.048	100mm	1.022.100	Ti.1.022.100
50mm	1.022.050	Ti.1.022.050			


ACL SCREW WITHOUT HEAD

8.0mm Dia.			9.0mm Dia.		
TOTAL LENGTH	CATALOGUE NO.		TOTAL LENGTH	CATALOGUE NO.	
	SS316L	TITANIUM		SS316L	TITANIUM
20mm	1.023.020	Ti.1.023.020	20mm	1.024.020	Ti.1.024.020
25mm	1.023.025	Ti.1.023.025	25mm	1.024.025	Ti.1.024.025
30mm	1.023.030	Ti.1.023.030	30mm	1.024.030	Ti.1.024.030
35mm	1.023.035	Ti.1.023.035	35mm	1.024.035	Ti.1.024.035
40mm	1.023.040	Ti.1.023.040	40mm	1.024.040	Ti.1.024.040
10.0mm Dia.					
TOTAL LENGTH	CATALOGUE NO.				
	SS316L	TITANIUM			
20mm	1.025.020	Ti.1.025.020			
25mm	1.025.025	Ti.1.025.025			
30mm	1.025.030	Ti.1.025.030			
35mm	1.025.035	Ti.1.025.035			
40mm	1.025.040	Ti.1.025.040			

ACL SCREW WITHOUT HEAD

8.0mm Dia.			9.0mm Dia.		
TOTAL LENGTH	CATALOGUE NO.		TOTAL LENGTH	CATALOGUE NO.	
	SS316L	TITANIUM		SS316L	TITANIUM
20mm	1.026.020	Ti.1.026.020	20mm	1.027.020	Ti.1.027.020
25mm	1.026.025	Ti.1.026.025	25mm	1.027.025	Ti.1.027.025
30mm	1.026.030	Ti.1.026.030	30mm	1.027.030	Ti.1.027.030
35mm	1.026.035	Ti.1.026.035	35mm	1.027.035	Ti.1.027.035
40mm	1.026.040	Ti.1.026.040	40mm	1.027.040	Ti.1.027.040
10.0mm Dia.					
TOTAL LENGTH	CATALOGUE NO.				
	SS316L	TITANIUM			
20mm	1.028.020	Ti.1.028.020			
25mm	1.028.025	Ti.1.028.025			
30mm	1.028.030	Ti.1.028.030			
35mm	1.028.035	Ti.1.028.035			
40mm	1.028.040	Ti.1.028.040			

DHS / DCS SCREW



TOTAL LENGTH	CATALOGUE NO.		TOTAL LENGTH	CATALOGUE NO.	
	SS316L	TITANIUM		SS316L	TITANIUM
50mm	1.029.050	Ti.1.029.050	100mm	1.029.100	Ti.1.029.100
55mm	1.029.055	Ti.1.029.055	105mm	1.029.105	Ti.1.029.105
60mm	1.029.060	Ti.1.029.060	110mm	1.029.110	Ti.1.029.110
65mm	1.029.065	Ti.1.029.065	115mm	1.029.115	Ti.1.029.115
70mm	1.029.070	Ti.1.029.070	120mm	1.029.120	Ti.1.029.120
75mm	1.029.075	Ti.1.029.075	125mm	1.029.125	Ti.1.029.125
80mm	1.029.080	Ti.1.029.080	130mm	1.029.130	Ti.1.029.130
85mm	1.029.085	Ti.1.029.085	135mm	1.029.135	Ti.1.029.135
90mm	1.029.090	Ti.1.029.090	140mm	1.029.140	Ti.1.029.140
95mm	1.029.095	Ti.1.029.095	145mm	1.029.145	Ti.1.029.145

DHS / DCS COMPRESSION SCREW



THREAD DIA.	CORE DIA.	HEAD DIA	HEXAGONAL SOCKET	TPI
3.80	2.80	9.0	3.5	32.0
TOTAL LENGTH	CATALOGUE NO.			
	SS316L	TITANIUM		
36mm	1.030.036	Ti.1.030.036		

4.50mm MALLEOLAR SCREW



THREAD DIA.	CORE DIA.	HEAD DIA	HEXAGONAL SOCKET	TPI
4.50mm	3.10mm	8.0mm	3.50mm	14
TOTAL LENGTH	CATALOGUE NO.		TOTAL LENGTH	
	SS316L	TITANIUM		
25mm	1.031.025	Ti.1.031.025	12.0mm	
30mm	1.031.030	Ti.1.031.030	15.0mm	
35mm	1.031.035	Ti.1.031.035	18.0mm	
40mm	1.031.040	Ti.1.031.040	20.0mm	
45mm	1.031.045	Ti.1.031.045	22.0mm	
50mm	1.031.050	Ti.1.031.050	24.0mm	
55mm	1.031.055	Ti.1.031.055	26.0mm	
60mm	1.031.060	Ti.1.031.060	28.0mm	
65mm	1.031.065	Ti.1.031.065	30.0mm	
70mm	1.031.070	Ti.1.031.070	32.0mm	
75mm	1.031.075	Ti.1.031.075	34.0mm	
80mm	1.031.080	Ti.1.031.080	36.0mm	

HERBERT SCREW



THREAD DIA. "A"	THREAD DIA.	SHAFT DIA.	HEXAGONAL SOCKET	CANNULATION
4.50mm	3.50mm	2.5	2.50	1.3
TOTAL LENGTH	CATALOGUE NO.		THREAD LENGTH A/B	
	SS316L	TITANIUM		
14mm	1.032.014	Ti.1.032.014	5.30/5.30mm	
16mm	1.032.016	Ti.1.032.016	5.30/5.30mm	
18mm	1.032.018	Ti.1.032.018	6.30/6.30mm	
20mm	1.032.020	Ti.1.032.020	6.30/6.30mm	
22mm	1.032.022	Ti.1.032.022	6.30/6.30mm	
24mm	1.031.024	Ti.1.032.024	6.50/6.30mm	
26mm	1.031.026	Ti.1.032.026	7.0/6.50mm	
28mm	1.031.028	Ti.1.032.028	7.0/6.50mm	
30mm	1.031.030	Ti.1.032.030	7.0/6.50mm	